


Singfiller[®]

— ENDLESS PURSUIT OF BEAUTY —



Singfiller® Cross-linked Sodium Hyaluronate Gel for Plastic Surgery




Certification: CE, NMPA, CFS, ISO13485, ISO9001

#Top HA Filler for Nasolabial Folds Reduction

- 1 Free from animal proteins
- 2 Fermentation technology
- 3 BD pre-fillable syringe

Application: Used for plastic surgery to remove wrinkles or Volume filling:



Singfiller® Production Process



Raw materials



Filling



Sterilization



Labeling



Packaging



Delivery

Where to use



Anti wrinkle



Lip filling




Nose enhancement



Breast enlargement

Principle of crosslinking



Note:

- ① Lower temperature by adding cross-linking agent, increase the stirring force, so that cross-linking agent and sodium hyaluronate dispersion is relatively uniform
- ② Control the cross-linking degree of crosslinked hyaluronic acid by controlling the reaction amount of hyaluronic acid and crosslinking agent
- ③ Access to products with cross-linked homogeneity through the rheological detection



Package



Prefilled BD
syringe



Vial



PACKAGE

In Blister



OEM





Global Network





Exhibitions



ARAB HEALTH

2019.1.28 – 31

2019.1.31 – 2.2
IMCAS
International Master Course on Aging Science

Dubai Derma
Skin Health is our Concern

2019.3.18 – 20

AMWC
17th AESTHETIC & ANTI-AGING
MEDICINE WORLD CONGRESS

2019.4.4 – 6

AFRICA HEALTH

2019.5.28 – 30

EXPO+MED 2019
JUNE 5th TO 7th
WTC • MEXICO CITY

2019.6.5 – 7

MEDICAL FAIR Bangkok Thailand
THAILAND 2019.9.11 – 13

SINGCLEAN
www.singclean.net
Exhibition

中国进出口商品交易会
CHINA IMPORT AND EXPORT FAIR
Since 1957

2019.10.31 – 11.4


MEDICA

2019.11.18 – 21
WORLD FORUM FOR MEDICINE
DÜSSELDORF, GERMANY
18–21 NOVEMBER 2019

Member of **MEDICAlliance**




Advantages







Parameter comparison





Viscoelasticity

—strong cohesion, perfect shaping







good viscosity—strong cohesion, hard to shift and spread
good elasticity—perfect shaping, deformation resistance



Particle size


——well-distributed particles, better shaping performance





Degree of swelling

——suitable swelling degree





- too high → easier to cause swelling;
- too low → easy displacement, causing foreign body sensation.



Push force

——uniform, stable, easier to inject






Feedback from Doctor



Effectiveness



Singfiller Indication



Type	Fine		Derm		Deep		Deeper		Sub-Q
Specification	1ml	/	1ml	2ml	1ml	2ml	1ml	2ml	2ml
Singfiller									
Indication	<p> Fine Derm Wrinkles across forehead Deep Deeper Temple Fine Frown lines Fine Derm Crow's feet Derm Deep Earlobe Fine Derm Nose Deep Deeper Cheek Deep Deeper Nasolabial folds Fine Derm Corners of mouth Fine Lip volume Derm Deep Chin filling Fine Derm Neck Lines Sub-Q Breast & Buttock </p>								



Injection Before & After



Before injection



1 month after
injection



6 months
after injection

SINGCLEAN



THANKS

